

# GSR2™

## **Biofeedback Relaxation System**


---

**INSTRUCTION MANUAL**

---



# ATTENTION



To avoid damaging connection & voiding your warranty **NEVER** plug earphone/meter into black input jack.

LB0040-00

**Note:** All CD referenced in the manual for the GSR 2 Basic relaxation instructions and guidance can now be found at the following link in MP3 format.

[thoughttechnology.com/index.php/hardware/gsr-2.html](http://thoughttechnology.com/index.php/hardware/gsr-2.html)

# GSR 2<sup>®</sup>



**Thought Technology Ltd.**

5250 Ferrier, Suite 812,  
Montréal, Québec  
H4P 1L3 Canada

**REF**

**SA2000M**



- Type BF Equipment
- Internally powered equipment
- Continuous operation
- Follow instructions for use



### **WARNING**



- Do not use in the presence of a flammable anesthetic mixture with air or with Oxygen or Nitrous Oxide.

### **ATTENTION**



- Keep Dry. Not to be immersed in water.
- Do not expose to extreme weather conditions.



- This device may be susceptible to electromagnetic interference. If abnormal performance is observed, additional measures may be necessary to operate the device in a suitable environment
- Not for diagnostic purposes. Not defibrillator proof. Not for critical patient monitoring.



- Do not leave batteries in the device for long periods of time when the device is not being used.

## **INTENDED PURPOSE**

- Biofeedback & Relaxation

Manual #SA2035 Rev.9 (May 2019)

# Table of Contents

|  |           |
|--|-----------|
| <b>Our age of stress .....</b>                           | <b>5</b>  |
| <b>Body language .....</b>                               | <b>5</b>  |
| <b>Biofeedback, explained .....</b>                      | <b>6</b>  |
| <b>Listen to your skin .....</b>                         | <b>6</b>  |
| <b>Thoroughly modern GSR 2 .....</b>                     | <b>6</b>  |
| <b>How the GSR 2 works .....</b>                         | <b>7</b>  |
| <b>Three Simple Relaxation Exercises .....</b>           | <b>10</b> |
| <b>Two advanced visualization exercises .....</b>        | <b>14</b> |
| <b>Benefits of the GSR 2 in your everyday life .....</b> | <b>14</b> |
| <b>Upkeep of the GSR 2 .....</b>                         | <b>15</b> |
| <b>Battery .....</b>                                     | <b>15</b> |
| <b>Optional Equipment .....</b>                          | <b>16</b> |
| <b>Specifications .....</b>                              | <b>19</b> |
| <b>Product Accessories .....</b>                         | <b>20</b> |
| <b>Warranty .....</b>                                    | <b>21</b> |
| <b>Returning Defective Equipment .....</b>               | <b>22</b> |
| <b>Placing Orders .....</b>                              | <b>22</b> |

## **Our age of stress**

It is ironic that today, when technology and science have provided our society with so much, when so many of us eat better, live better, enjoy more leisure and comfort than any other people in history, we suffer from such widespread stress and tension.

Unless we live protected lives as virtual hermits on some island paradise, stress is very much with us - unfortunately for many, far too much.

The dangers of excessive stress to our mental and physical health are well documented, the all-too-frequent answer being a reliance on alcohol or tranquilizing drugs. Fortunately, people are seeking new ways to cope, and are making everything from active sports to personal development programs, meditation and biofeedback training part of their everyday routines.

The ability to calm down and relax deeply is a powerful antidote to the effects of stress. The problem is that although all of us are born with this ability, it's something we often have to re-learn.

## **Body language**

When you're under extreme stress, your body copes in ways that are very noticeable. Your heart pounds, your hands perspire, muscles tense up, extra adrenaline is released.

But not all stress is so extreme. Unfortunately, what really wears us down are lower levels of stress that occur constantly and that most of us accept as inevitable. Our bodies respond to this type of stress with much lower signals that mostly go unheeded. Now, biofeedback is enabling us to observe and modify these signals.

## **Biofeedback, explained**

Biofeedback is an exciting field that combines physiology and psychology. The biofeedback technique uses electronics to detect and amplify internal body activities too subtle for normal awareness.

Once you are able to observe these activities, you can then try to modify them. Because you learn immediately whether you have succeeded or not, you can gradually teach yourself to produce the desired results.

## **Listen to your skin**

One of the principal ways your body reacts to tension and stress is through your skin and what scientists call galvanic skin resistance.

This galvanic skin resistance, or GSR, is a reflection of variations in your sweat gland activity and pore size, both of which are controlled by the sympathetic nervous system. When you become excited, frightened or disturbed to any degree, the system activates chemical and physical changes all through your body. Your GSR level also changes.

Your skin resistance increases when you are calm and relaxed. When you tense up even slightly, your skin resistance decreases.

By letting you hear those increases and decreases, the GSR 2 monitor can give you a new awareness of the internal state of your mind and body, and help you develop the ability to control tension and stress.

## **Thoroughly modern GSR 2**

Psychologists have been monitoring the electrical resistance of the skin to gauge nervous system activity for over 65 years. The first instruments used featured an early-model electric meter called a galvanometer, after the 18th



century Italian scientist, Luigi Galvani. He also gave his name to GSR.

By the 1950's, GSR machines were vastly improved, though still cumbersome and complicated by today's standards. They were available almost exclusively to medical practitioners and research personnel.

Next came the GSR 2's immediate ancestors: biofeedback instruments that were light, battery-powered and solid state, but still cumbersome. The GSR 2 is a great step forward.

Created and developed by clinical psychologists and engineers working in the field of stress reduction, the GSR 2 has done away with finger attachments. Its electrodes are an integral part of the case.

The GSR 2 is automatic. No switch is needed to turn it on or off. It can be used just about anywhere.

Most important, the GSR 2 is so reliable, durable, affordable and simple to use, it has brought the benefits of biofeedback training within reach of virtually everyone.

Now in use around the world, the GSR 2 has been well accepted by leading medical, educational, athletic and social institutions.

## **How the GSR 2 works**

To turn the unit on, place your fingers on the sensing plates. Turn the dial *very slowly* until you hear a tone. When you start to become calm and relaxed, your skin resistance increases and the tone lowers. When you become upset or agitated, your skin resistance decreases and the tone goes higher.

Put very simply, you learn to relax by learning how to diminish the tone of the GSR 2 to the lowest point possible.

Read the instructions thoroughly and then you're ready to enjoy your first skin biofeedback experience.

The GSR 2 will shut off automatically as soon as you remove your fingers from the sensing plates.

## **Instructions**

To get the most from your GSR 2, read the instructions before use. You will also find it helpful to listen to your CD a few times.

1. Try to set aside the same relaxation time for yourself every day. How long you spend is up to you. Start with 10-20 minutes, or as long as you feel comfortable.
2. Find a quiet spot where you won't be interrupted or distracted. Average room temperature is best for operating the unit. Relax by loosening any tight clothing and removing your shoes.
3. Sit in a comfortable position at a desk, in an armchair or anywhere there is arm and elbow support. Lying down is not recommended unless your purpose is to fall asleep.
4. The GSR 2 can be used with either hand. However, using your left hand frees your thumb for easier dialing. To find the right position for your arm, let it rest naturally on its side, as if you were holding a glass. The unit should rest easily in your hand.
5. Place your first two fingers limply on the sensing plates; your thumb goes underneath the unit. Never press down or squeeze. Only your natural touch is needed, as the unit is designed to allow full and even contact, even if your fingers do not cover the complete length of the plates.

6. If your hand needs steadying, flip the strap up. The strap is not necessary; however, it is recommended that you use it. Use it to stop hand movement and help you maintain a light, even pressure on the plates. To adjust band tightness, gently pull the band out on both sides. Many people prefer to allow the GSR 2 monitor to rest upside down in the palm of their hand. (See illustration 6a). Some people use the remote electrodes, described on page 18 of this manual, to permit stable monitoring even while moving their hand.
7. With your fingers in position, turn the dial all the way towards your wrist with the thumb of your other hand. Slowly turn it back up until you hear a pleasant, low tone. If you turn the dial too far up, the tone will become too high-pitched for you to hear. Should that happen, simply turn the dial back slowly. Wait about 20 seconds for the tone to stabilize, then check once more to make sure you're at the low end of the tonal range. Once your tone level is established, leave it there.
8. Now you're ready to begin. With your hand in place and the tone set low, sit quietly with your eyes closed. Relax your muscles and breathe evenly. Let go of all distracting thoughts and allow your mind to clear. As you become more relaxed, the tone will gradually decrease. Don't expect instant results. Trying too hard can actually raise the tone.

There is usually about 2-5 seconds delay before your pores react, so don't expect instant tonal changes. Also, bear in mind that because you relax more slowly than you tense up, the tone will take longer to fall than to rise.

As you become more skillful at relaxing, the tone may diminish until you can't hear it. Just reset the

dial slightly higher and continue as before. If the tone ever goes too high to hear, simply reset the dial to a lower starting point. The numbers on the dial are marked 1... to 9... You may like to note your setting before and after use.

With the GSR 2 tone at the same pitch, a higher number indicates relaxation. For example, if following your biofeedback session, the dial reads 2, yet later when you reset the dial to produce the same pitch, it reads 3... it indicates a decrease in tension. It is meaningless to compare your setting with other people since each person's skin is as unique as his/her fingertips. Remember, your starting and finishing point is of little importance, rather it is the feeling of physical and mental relaxation which is the goal.

*People who learn biofeedback control cannot always say how it is done. But be assured that before long, if you practice regularly, you will work out your own method. Don't be impatient. You have an important goal in mind: not only to bring the tone down, at will, as low as possible, but also to teach yourself to cope with stress through relaxation.*

## **Three Simple Relaxation Exercises**

*The following exercises demonstrate how the GSR 2 instantly and continuously monitors your progress as your mental and physical states go from tense to relaxed.*

*NOTE: Throughout these exercises, try to maintain even, relaxed breathing. (This is well described in your recording.)*

## 1) Basic Tension Relief

The arrows indicate three major stress points in you body. Set your GSR 2 at a low tone and try to encourage relaxation in these areas by doing the following:

- Let your shoulders fall.
- Let your jaw drop. As it does, your mouth will open slightly.
- Relax the muscles around your eyes and forehead.

## 2) Muscle Relaxation

\* Instructions by Edmund Jacobson MD

Prepare for a biofeedback session with the GSR 2. Place it in a position that you can easily access, say near the right hand. Sit or lie down to use.

### Before using the GSR 2:

- Place both hands gently down on a flat surface, palms down.
- Place the fingers of your right hand on top of the soft fleshy muscle mass just after the elbow bend on the top of the forearm.
- Make a tight fist with the left hand. Notice the muscles under your right fingers tighten as you hold a tight fist. This is you making an EFFORT to hold a tight fist.
- Loosen the tight left fist and notice the tight muscles under your right fingers go limp at the same time. That is you stopping the EFFORT of tightening the fist. Be sure that you are using your fist and not the muscles in your neck and shoulder.
- Return both hands to the palm down starting position.

## Using the GSR 2:

- Place your right hand on the GSR 2 as if it were relaxed with your fingers under the finger band. Adjust the tone to a pleasant low pitch.
- Again, make a tight fist with your **left** hand for several seconds, just as you previously did. Be sure to breathe evenly and deeply during the entire session.
- This time focus on the sensation in the left forearm, about at the same place that you had placed your **right** fingers to feel the muscles tense. There you will find the **sensation of you tensing** that we want you to learn to recognize everywhere in your body.
- Let go of the tight left fist. Notice the **sensation of you tensing** goes away.
- Again, make a tight fist with the **left** hand. This time notice what happens to the tone of the GSR 2. If the tone increases, it indicates you have increased SNS arousal.
- Gradually ease the tension in your **left** fist and arm. Notice the easing of the **sensation of tension** in the **left** forearm, and anywhere else in your body where you tensed in concert while tightening the fist. Notice, that in general, the tone of the GSR 2 will gradually decrease to a lower range as you yourself decrease your tension and experience a greater sense of relaxation.
- Relax for a minute or two, and just observe the GSR2's tone.
- Again, tighten the left fist, but this time only about one half as firmly as you did the first time. Again, observe any increases in your GSR2's tone. Notice **the sensation of tenseness** has returned to the same place just below your left elbow bend, but the sensation should be less distinct or intense.

- Try to let the tightness in your fist and fingers go bit by bit, or step by step and see if you can notice the **sensation of tenseness** fading step by step.
- See if you can reach both a lower **sensation of tenseness and GSR tone** than when you started.

If you can learn to perceive your tension by recognizing the **sensation of tension**, then you are starting on your goal to learn progressive relaxation. To generalize your new ability to relax to other parts of the body, try to find the **sensation of tenseness**, in other parts of your body and let go in the same stepwise manner.

### 3) Breathing Relaxation

- Keep your hand relaxed on the GSR 2 and wait until the tone has stabilized in the low range.
- Relax quietly with your eyes closed and focus on your breathing. Each time you inhale, fill your diaphragm by allowing your stomach to expand and letting this expansion roll up to your chest. (Refer to your recorded relaxation exercise if necessary.) Think of the tone as a clock ticking in the background as you focus your attention on your breathing.
- Exhale slowly and completely, each time thinking the word “calm”, or any other word you would prefer to use. Stretch out the word so it becomes CAAAALMMMM.
- If thoughts arise or your attention wanders, simply focus on your breath when you realize what has happened.
- This exercise may help you learn to let go of distracting or tension-producing thought patterns.

## **Two advanced visualization exercises**

1. Visualization is a useful relaxation technique. When you have mastered the previous exercises, raise the tone of the GSR 2 to a higher starting level and teach yourself to lower it to your previous level. Visualize and experience in detail a pleasant situation. Observe its effects on the tone. Then visualize an unpleasant or stressful situation. If this raises the tone, concentrate on bringing it down again by switching back to the pleasant scene.
2. This one is more difficult, but very rewarding if you master it. After you have succeeded in lowering the tone, visualize an anxiety-producing situation and try to maintain the low tone by countering your anxiety with relaxation. The purpose of this exercise is to encourage you to relax when confronted with a real, stressful situation.

## **Benefits of the GSR 2 in your everyday life**

1. Once you are able to put yourself into a relaxed state, you have a valuable skill that you can apply to your everyday life. Your increased sensitivity to yourself will sharpen your awareness of how you respond to stressful situations, allowing you to remain calm, even without your unit.
2. Deep relaxation gives the mind and body a chance to recuperate. When your body is released from stress, you can enjoy increased energy.
3. The GSR 2 is often used in conjunction with personal development programs such as autogenic training, yoga and meditation. It acts as an internal guide, providing feedback on your state of focus and detachment, as well as on your improving skills.
4. GSR biofeedback skills can be put to good use by sports participants. Athletes do better in tense



situations if they can channel their excess nervous energy into performance and allow their natural abilities to surface unimpeded.\*

5. Likewise, learning GSR biofeedback skills can be of great benefit to students and people working in high-pressure jobs.

\* Mind Over Muscle™ training program tapes and workbook for athletes (T2861).

## **Upkeep of the GSR 2**

The GSR 2 is made of high quality plastic and polished stainless steel. It is extremely durable, with no special care needed. However, as with any electronic instrument, it should not be handled abusively. To clean, wipe with a damp cloth and a little mild soap, if needed. **DO NOT IMMERS**E IN WATER. The finger plates should merely be wiped with a soft, clean cloth.

## **Battery**

The GSR 2 has been designed to minimize the use of battery current. The unit should function for the shelf life of the battery (6 months to 2 years). A decrease in volume and sensitivity indicates that the battery needs changing.

### **To change:**

Remove the two screws and carefully separate the unit into halves. Pry off the battery snaps, then remove the battery. Replace it with a high-quality 9v cell, preferable alkaline. Carefully fit the unit back together. When closing the unit, turn the screws only until resistance is felt.

# Optional Equipment

## 1. Reverse Temperature Sensor

The reverse temperature sensor adds a completely different function to your GSR 2 and lets you experience a fascinating new area of biofeedback technology - temperature biofeedback. Generally, increasing levels of stress lessen the blood flow to your hands and feet. This causes decreases in hand and foot temperatures. With the monitoring help of the thermistor, you can learn to control that temperature. It's amazing when you realize what you're basically controlling is the dilation and constriction of blood vessels in your arms and legs.

The temperature sensor can also be used with the relaxation exercises described earlier in this manual.

### **To use:**

1. Plug the temperature sensor into pin-protected input jack of the GSR 2.
2. Place the tip of the sensor on the pad of a finger (either hand) or on any toe. Use a different location each time you train to promote general hand or foot warming.
3. Wrap the Velcro® around your finger or toe to hold the sensor in place, using only enough pressure to hold it there. Do not apply pressure or inaccurate readings may result.
4. *Important:* Do not touch GSR 2 sensing plates while using the temperature sensor, as this will give you a false reading. Remember, the idea of relaxation with the reverse temperature sensor is to lower the tone, which indicates the hands are warming.
5. Set the dial of the GSR 2 to a pleasant low tone, or center the meter if you're using it. To lower the tone,

relax completely and think of your hand and arm, or foot and leg, as heavy and warm.

6. Note: The temperature sensor must be unplugged when not in use to avoid battery drain.



You may use a home thermometer to show a rise or decline of hand or foot temperature.

## **Visualization Exercise for Temperature Biofeedback**

1. Plug the temperature sensor into the pin-protected jack, wrap the Velcro® around your finger or toe, and get comfortable.
2. Close your eyes.
3. Focus on your breathing. Breathe evenly and deeply for several minutes.
4. Now visualize a warm place where you are relaxed and comfortable - on the beach, in a sauna or in the hot sun.
5. When the scene is fixed in your mind, think of your warm hands and feet. Feel the sensation and passively observe it.
6. When you begin to relax, your hands and feet will get warmer and the tone of the GSR 2 will fall. Passively observe this as well.
7. If you are successful at raising your temperature, pay attention to the internal sensations that accompany

the increase. In that way, you may learn to achieve more direct control over hand and foot temperature variations without any electronic biofeedback help.

## 2. Dual Sensitivity Meter

Plugging the meter into the GSR 2's 2.5 mm jack automatically turns off the tone and provides visual feedback on either temperature or skin resistance. The meter faces are marked with division of  $.1^{\circ}\text{F}$  and  $.05^{\circ}\text{C}$ . The dual sensitivity meter switched to the "1" position provides a maximum reading of  $\pm 2^{\circ}\text{F}$  and  $\pm 1^{\circ}\text{C}$ . The "1/2" position doubles this range to  $\pm 4^{\circ}\text{F}$  and  $\pm 2^{\circ}\text{C}$ .



When monitoring with the reverse temperature sensor, the needle moves counter-clockwise with warming.

For monitoring skin resistance, place your fingers on the sensing plates, or attach the remote electrodes, and slowly turn the dial on the GSR 2 until the meter needle centers. The meter moves counterclockwise as you relax. The numbers have no specific meaning in this mode, but can be used to note changes.

*For users whose GSR changes rapidly, switching to the 1/2 sensitivity setting on the dual sensitivity meter doubles the training range and decreases frequent resetting of the dial.*

### 3. Electrodes

The electrodes plug into the pin-protected input jack of the GSR 2. They can be used if you have difficulty holding your fingers still on the GSR 2's fingerplates. (Cat. No. SA2050)



#### Specifications

Skin resistance range 1,000 ohms - 3,000,000 ohms.  
Variable frequency range 0 to 40,000 Hz.  
9 volt battery. Battery current in use < 2.0 mA.  
2.5mm output jack for 8 ohm earphone and output meter.  
pin-protected input jack for thermistor and electrodes.

Operating conditions:

Temperature: +5°C to +40°C

Atmospheric Pressure: 700hPa to 1060hPa

Relative Humidity (non-condensing): 15% to 93%

Storage and transportation conditions:

Temperature: -25°C to +70°C

Relative Humidity (non-condensing): 93%

# Disposal



Appropriate disposal of the device and sensors should be done in accordance with accepted medical practice and any applicable local, state and federal laws and regulations.

## Product Accessories

### **T2001M .....GSR Unit**

- SA2000M .....GSR 2 Unit
- SA3409 .....MyoTrac/Dac Earphone 8 Ohms
- SA2032 .....CD GSR2 English
- SA2035.....Manual GSR2 English

### **T2120M ..... GSR 2/Temp 2X**

- SA2000M .....GSR 2 Unit
- SA3409 .....MyoTrac/Dac Earphone 8 Ohms
- SA2032 ..... CD GSR2 English
- SA2035 .....Manual GSR2 English
- SA2050M .....Remote Electrodes
- SA2076 .....Dual Sensitivity Meter
- SA2090M .....Reverse Temperature Sensor
- SA2007 .....Thermometer (Not for European sale)

## Warranty

The GSR 2 and all equipment, including optional items, are guaranteed to be free from defects in material and workmanship for one year from the date of purchase. In the unlikely event that repair is necessary, Thought Technology Ltd. will repair or replace your GSR 2 free of charge. Return the unit postage prepaid and insured, with proof of purchase.

This warranty does not apply to damage, whether incurred through accident, alteration or abuse.

No other warranty is expressed or implied.

*The GSR 2 is battery-operated. Batteries are not covered under this warranty. Remove dead batteries promptly to prevent corrosion damage.*

## Returning Defective Equipment

1. Send the unit(s) **postage prepaid** and **insured**, with proof of purchase to one of the addresses below.
2. If you are shipping from anywhere outside of Canada or from the USA to Canada, mark the package "**Goods to be repaired - Made in Canada**" to avoid unnecessary customs charges.  
Overseas repairs must write Thought Technology broker name and account number **Livingston International – 133461** outside the package and in the commercial invoice.
3. All customs and duty charges will be billed to you if incurred by sending the unit to the **wrong** address.
4. Provide a **detailed** description of the problem you are experiencing, and your telephone/fax number.

### Ship insured to:

|  |   |
|--|---|
| In the USA:<br><b>Thought Technology Ltd.</b><br><b>Cimetra LLC</b><br><b>8396 State Route 9</b><br><b>West Chazy, New York</b><br><b>12992, USA</b> | In Canada and other countries:<br><b>Thought Technology Ltd.</b><br><b>5250 Ferrier, Suite 812,</b><br><b>Montréal, Québec</b><br><b>H4P 1L3 Canada</b> |
|--|---|

[www.thoughttechnology.com](http://www.thoughttechnology.com)

## Placing Orders

To order call: **1 (514) 489-8251**

FAX orders: **1 (514) 489-8255**



|                                  |       |
|----------------------------------|-------|
| Index                            |       |
| band .....                       | 9     |
| Basic Tension Relief ....        | 11    |
| battery .....                    | 15    |
| Biofeedback.....                 | 6     |
| biofeedback instruments          | 7     |
| Breathing Relaxation ...         | 13    |
| cleaning.....                    | 15    |
| Deep relaxation .....            | 14    |
| Defective Equipment....          | 22    |
| dial .....                       | 7     |
| Dual Sensitivity Meter ..        | 18    |
| false reading.....               | 16    |
| finger band .....                | 12    |
| galvanic skin resistance .       | 6     |
| galvanometer.....                | 6     |
| GSR .....                        | 6     |
| GSR 2 monitor .....              | 9     |
| Muscle Relaxation.....           | 11    |
| nervous system.....              | 6     |
| Orders.....                      | 22    |
| Relaxation Exercises....         | 10    |
| sensing plates .....             | 7     |
| skin resistance .....            | 6, 18 |
| stress.....                      | 5     |
| temperature biofeedback<br>..... | 16    |
| temperature sensor .....         | 16    |
| tension .....                    | 6     |
| tonal changes .....              | 9     |
| tone .....                       | 7     |
| Visualization.....               | 14    |
| warranty.....                    | 21    |

The logo consists of two stylized, overlapping 'T' shapes. The top 'T' is white with a blue shadow, and the bottom 'T' is blue with a white shadow, creating a 3D effect.

# Thought Tech

**Fully guaranteed for one year.  
Sold in more than 70 countries around the  
world. Used in thousands of leading medical,  
educational and social institutions. The GSR  
2 Biofeedback Relaxation System is one of  
the fine products developed by:**

**Thought Technology Ltd.**

**5250 rue Ferrier, Suite 812**

**Montreal, Quebec, Canada H4P 1L3**

**U.S. Pat. Nos 3648686; D 38962**

**Canada Rd. 1975 Thought Technology Ltd. Foreign  
Design Reg'ns: Australia, France, Great Britain,  
Switzerland.**

**Other design registrations and patents pending.**



1975, 1977, 1978

International Design Award

Geneva, Switzerland. Nuremberg, Germany

