



Thought Technology Ltd.

8205 Montreal-Toronto Blvd. Suite 223, Montreal West, Quebec, H4X 1N1 Canada
 Tel: 800-361-3651 (514) 489-8251 Fax: (514) 489-8255
 Email: mail@thoughttechnology.com
www.thoughttechnology.com

Recommended Computer Requirements for BioGraph Infiniti System:

Processor	Quad-core CPU, Intel i5
Operating systems (supported) (32 and 64 bit versions)	<ul style="list-style-type: none"> • Microsoft® Windows® 10, Pro, Enterprise • Microsoft® Windows® 8, Pro, Enterprise • Microsoft® Windows® 7 Home Premium, Professional, Enterprise or Ultimate with SP1 • Microsoft® Windows Vista® Home Basic, Home Premium, Business, Ultimate, or Enterprise with Service Pack 2
Memory (RAM)	3 GB (Vista® or Windows® 7/8/10)
Hard Disk space	50 GB or more, to save video and sessions data
Video / Graphics card	At least 1 GB of dedicated memory. (Example: ATI Radeon 5770, with 1280x1024 resolution and dual monitor support.)
Monitor	Two monitors which support resolution of 1280x1024 and higher.

Other necessary components:

- DVD-ROM drive.
- Adequate number of USB ports depending on the intended system usage (1-4 for Pro/Flex Infiniti encoders, 2 for webcams, 1 for TT AV Sync, mouse, keyboard, printer)
- Mouse or any other compatible pointing device.
- 32-bit compatible sound card & speakers.
- Compact Flash Reader (for use with compact flash card only).
- Up to two webcams with built-in microphone (if recording video). Frame rate 30fps, Microsoft® recommended.
- Microsoft® Office with Word® and Excel® (for report generation and printing).
- High-speed Internet access (for updating software and receiving technical support).

PLEASE NOTE:

- Tele-Infiniti CF (T9600) version 3.2.1 or higher is compatible with 64-bit Windows® operating systems. Version 3.1 or lower is compatible with 32-bit Windows®.
- Low speed CPUs (1.8 GB) may be sluggish when running screens with a large number of instruments per screen.
- If using 2 webcams for video recording, it is strongly suggested to use recommended PC requirements.
- If using TT-AV Sync functionality, please refer to possibly additional system requirements from appropriate suite manual.

The following apply to BioGraph Infiniti 6.1 and its application suites:

- BioGraph Infiniti is a 32-bit software program.
- BioGraph Infiniti software and application suites are installed for **All Users** on a per-machine basis.
- The BioGraph Infiniti system must be installed on the PC where the database is/will be located.
- BioGraph Infiniti is not designed for use in a networked environment