



BioGraph 6.2's exciting new features simplify your work so you can concentrate on teaching self-regulation skills to your clients: Infinitely expand your ability to engage and motivate your client with local or Web-based video content. Automatically and instantly fit your feedback screen to the available monitor size without fuss.

NEW FEATURES

VIDEO STREAM INSTRUMENT

Do biofeedback with video, stored on your PC or streamed from the Internet.

- Actions include shrink/enlarge, play/pause, blur, pixelate, volume loud/quiet and tunnel in/out.
- Custom select videos for each client from your library of available videos.
- When using videos stored on your PC, a bookmark is placed in the video at the end of a session, ready for the client's next visit.
- No editing required: Convert the DVD instrument on your favorite screen to a Video instrument.

AUTO-RESIZING

BioGraph 6.2 makes your life easier by giving you the option to automatically maximize feedback screens to fit your monitor, whatever its resolution.

SCREEN TEMPLATE EDITOR

Suite designers can move to the next level and create feedback screens that are channel set independent.

- You create a template by placing instruments on the screen and assigning a signal type to each instrument.
- Clinicians convert a template to any channel set and the editor automatically scans for appropriate channels to connect to each instrument.

