



The "SUITE" Deal Just Got Better!



Thought Technology is proud to launch two new specialized Suites: **EEG Suite and Physiology Suite**. The Suites take full advantage of functions and features of the powerful Biograph Infiniti software platform.

They provide complete control over your entire clinical process, including: assessment, training, reporting and demonstrating outcome. For each protocol, the Suites provide signal quality verification screens, multimedia feedback screens and quick report screens showing in-session trends & learning curves.

The EEG Suite contains 6 scripts and more than 100 screens (over 21 from expert clinicians) for 1 or 2-channel EEG on 1 or 2 monitors. It comes with both standard protocols and a variety of template screens for user-defined protocols:


- Standard Bands: Alpha, Alpha peak frequency, Alpha-Theta, Low-high Beta, Beta-Theta, SMR-Theta, Wide-band inhibit, etc.
- 3 User Bands: 1 reward (up), 2 inhibits (down)
- 6 User Bands: 2 rewards (up), 4 inhibits (down)
- Multiple Bands: 1 reward, 9 inhibits


The EEG Suite also offers many advanced tools for assessment and follow-up:

- Pre/Post-training baseline assessment scripts
- A 4-activity EEG assessment script: eyes-open, eyes-closed, and 2 tasks
- Cross-session trend tools to show your client's progress
- A set of hybrid EEG and peripheral integrated biofeedback screens

EEG Suite features include:

- On-the-fly cut-off frequency adjustments
- Impedance detection
- Auto-threshold instruments with user-specified target percentage of time
- Single and dual-monitor display screens
- Real-time epoch statistics for instant trend graphing
- Out-of-this world alpha-theta relaxation script





The Physiology Suite contains a complete set of standard biofeedback tools to get a total physiological perspective on your client's sympathetic and parasympathetic processes, including: Temperature, Skin Conductance, Surface EMG, Respiration and Heart rate variability (HRV).

The Physiology suite contains more than 80 screens and over 14 scripts, and provides you with advanced tools for assessment and follow-up:

- Pre/post baseline assessment
- SC verbal Relaxation exercise
- Training scripts (10 x 1 min trials)
- A 7 activity stress assessment script that takes your client through a sequence of stress/relax cycles and generates a report to rapidly identify hyper- or hypo-reactivity in each modality.

The two Suites offer 17 new animations, 15 MP3 (music & ambiance), 20 MIDI songs and many other new multimedia feedback options:

- DVD (4 action control: shrink, zoom, pause & sound)
- Morphing animation sequences, 3 stage puzzles
- Nature sound audio tracks
- MIDI splitter






Thought Technology Ltd.

2180 Belgrave Avenue Montreal, QC, H4A 2L8, Canada
Tel: 1-800-361-3651 • 514-489-8251 Fax: 514-489-8255
<http://www.thoughttechnology.com>

Physiology Suite features:

- Skin-conductance driven relaxation script with audio induction track
- User-defined pre-labelled event-markers (hot key activated)
- Real-time epoch statistics for instant trend graphing



MAR803-00

