

**European Community Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE).**

The information below identifies the different components and materials and the location of dangerous substances and preparations in the equipment. The information is posted on our website [www.thoughttechnology.com](http://www.thoughttechnology.com)

<b>Equipment</b>	<b>Individual Parts</b>	<b>Components &amp; Materials</b>	<b>Dangerous Substances &amp; Preparations</b>
GSR2	GSR2 unit Temperature sensor Remote electrode Dual meter Earphone	Lead Free Solder Printed circuit boards. SMT electronic components. Plastic enclosures. Foam packaging. Plastic insulation. Copper wiring. Steel electrodes. 9 Volt battery.	None
Myotrac	Myotrac unit Myoscan Sensor Earphone	Lead Free Solder Printed circuit boards. SMT electronic components. Plastic enclosures. Foam packaging. Plastic insulation. Copper wiring. Steel electrodes. 9 Volt battery.	None
Myotrac Infiniti	Myotrac Infiniti unit Headphone Power Adaptor Extender cables USB receiver	Lead Free Solder Printed circuit boards. SMT electronic components. Plastic enclosures. Cloth carry-case. Foam packaging. Plastic insulation. Copper wiring. Magnet Steel electrodes. 1.2 Volt batteries (4).	None
Procomp & Flexcomp Infiniti	Procomp & Flexcomp units Headphone Extender cables HRBVP sensor Skin Conductance sensor Respiration sensor Temperature sensor EEG sensor EKG sensor Myoscan sensor USB receiver	Lead Free Solder Printed circuit boards. SMT electronic components. Plastic enclosures. Cloth carry-case. Foam packaging. Plastic insulation. Copper wiring. Steel electrodes. 1.5 Volt batteries (4).	None

Equipment	Individual Parts	Components & Materials	Dangerous Substances & Preparations
Procomp2	Procomp2 unit Headphone Extender cables HRBVP sensor Skin Conductance sensor Respiration sensor Temperature sensor EEG sensor EKG sensor Myoscan sensor USB receiver	Lead Free Solder Printed circuit boards. SMT electronic components. Plastic enclosures. Foam packaging. Plastic insulation. Copper wiring. 1.5 Volt battery.	None
U-Control	U-Control unit Earphone Extender cables	Lead Free Solder Printed circuit boards. SMT electronic components. Plastic enclosures. Plastic carry-case. Plastic insulation. Copper wiring. 9 Volt battery.	None
eVu TPS	eVu TPS Sensor Power Adapter TPS Strap	Lead Free Solder Printed circuit board SMT electronic components Plastic enclosures Steel electrodes Plastic label 3.7V Lithium Polymer Battery knit polyester fabric	None
MyOnyx	MyOnyx device Pressure Sensor Vaginal Pressure Probe Patient Drive Cable EMG/STIM Cables Electrodes Din to Snap Cable Adapter Power Adapter	Lead Free Solder Printed circuit boards SMT electronic components Plastic enclosures 3.7V Lithium Polymer Battery Aluminum ring Copper wiring Silicon Probe Steel electrodes	None