



Thought Technology Ltd.

8205 Montreal-Toronto Blvd. Suite 223, Montreal West, Quebec, H4X 1N1 Canada
 Tel: 800-361-3651 (514) 489-8251 Fax: (514) 489-8255
 Email: mail@thoughttechnology.com
www.thoughttechnology.com

Recommended Computer Requirements for BioGraph Infiniti System:

Processor	Quad-core Intel i5 or a quad-core equivalent
Operating systems (supported) (32-and 64-bit versions)	<ul style="list-style-type: none"> • Microsoft® Windows® 10 Home, Pro, Enterprise • Microsoft® Windows® 8, Pro, Enterprise • Microsoft® Windows® 7 Home Premium, Professional, Enterprise, or Ultimate with SP1
Memory (RAM)	Minimum 6 GB (Windows® 7/8/10 64-bit) or 4 GB (Windows® 7/8/10 32-bit)
Hard Disk space (Primary C drive)	50 GB or more, to save video and sessions data on the primary drive (recommended to be no less than 200 GB)
Video / Graphics card	At least 1 GB of dedicated memory with 1280x1024 resolution and higher for 32-bit systems. Not necessary for recent generations of Intel 64-bit systems with a minimum of 8 GB RAM.
Monitor	Resolution of 1280x1024 and higher; option for two monitors

Other necessary components:

- Adequate number of USB ports depending on the intended system usage (1-4 for ProComp or FlexComp Infiniti encoders and accessories)
- 32-bit compatible sound card & speakers.
- Microsoft® Office with Word® and Excel® (for report generation and printing).
- High-speed Internet access (for streaming videos, updating software, and receiving technical support)

Optional components, depending on intended system usage:

- DVD-ROM drive for playing movies from a DVD; note: many non-DVD video options are available.
- Compact Flash Reader (for use with compact flash card only).
- Up to two webcams with built-in microphone (if recording video). Frame rate 30fps, Microsoft® recommended.
- Bluetooth receiver/radio, if using TPS sensor and other Bluetooth-enabled devices.

PLEASE NOTE:

- Tele-Infiniti CF (T9600) version 3.2.1 or higher is compatible with 64-bit Windows® operating systems. Version 3.1 or lower is compatible with 32-bit Windows®.
- Low speed CPUs (1.8 GHz) may be sluggish when running screens with a large number of instruments per screen.
- If using 2 webcams for video recording, it is strongly suggested to use recommended PC requirements.
- If using TT-AV Sync functionality, please refer to possibly additional system requirements from appropriate suite manual.

The following apply to BioGraph Infiniti 6.6 and its application suites:

- *BioGraph Infiniti is a 32-bit software program. HOWEVER, we recommend 64-bit version of Windows OS to be able to have 4+ GB of RAM.*
- *BioGraph Infiniti software and application suites are installed for **All Users** on a per system basis.*
- *The BioGraph Infiniti system must be installed on the PC where the database is/will be located.*
- *BioGraph Infiniti is not designed for use in a networked environment.*